**IIW's Analytical instruments** 

# NEVINEES.

Thread-forming properties (spinnability) measurement instrument.

## Simple & Easy

Measuring stretching properties of liquids







IMI-0901 IMI-0501



## Spinnability measurement Instrument NEVA METER

## **Easy measurements of spinnability**

## ■ Features

- · Simple alternative instrument for measuring liquid spinnability.
- · Requires a very little amount of sample solution.
- · No need of any reagent.
- · Repeats measurements (clotting test, etc.)
- · Excellent for urgent testing needs.
- · Minimal damage to sample/specimen.
- · No running costs. Optimal for screening tests.
- On-board electronic scale measures tensile load with an accuracy of 0.1 mg.(HI0901)

## ■ Specifications (overview)

| Products         |                        | IMI-0901   | IMI-0501   |
|------------------|------------------------|--|--|
|                  |                        |  | NA PARAMETER AND A PARAMETER A |
| Electronic scale |                        | Weighing capacity 220g   | _  |
|                  |                        | Data transfer 92bps  | _  |
|                  |                        | Accuracy 0.1mg   | _  |
| Computer         |                        | Compact factory computer   | _  |
|                  |                        | 19" LCD  | _  |
|                  |                        | External 320GB HDD   | _  |
| NEVA®METER       | Measurable length      | Effective stroke: 100.0mm  | Effective stroke: 100.0mm  |
|                  | Measuring speed        | Standard speed: 10.0mm/sec ±1%   | Standard speed: 10.0mm/sec ±1%   |
|                  |                        | Variable range: 4.0 ∼ 100.0mm/sec  | Variable range: 4.0 ∼ 100.0mm/sec  |
|                  | Measurement accuracy   | Minimum display value: 0.01mm  | Minimum display value: 0.01mm  |
|                  |                        | Measurement accuracy: ±0.05mm  | Measurement accuracy: ±0.05mm  |
|                  |                        | Accuracy of repeated measurements with standard liquid: ±0.02mm  | Accuracy of repeated measurements with standard liquid: ±0.02mm  |
|                  | Specimen required      | $\geq$ 0.06m $\ell$  | $\geq$ 0.06m $\ell$  |
|                  | Probe and testing dish | Standard probe (3-mm diameter, flat)   | Standard probe (3-mm diameter, flat)   |
|                  |                        | Standard testing dish (0.08 ml drop)   | Standard testing dish (0.08 ml drop)   |
|                  | Cleaning solution      | Ethanol for disinfection   | Ethanol for disinfection   |
|                  | Display                | Liquid crystal display (16 characters × 2 lines, interactive)  | Liquid crystal display (16 characters × 2 lines, interactive)  |
|                  | Power                  | AC100V 50/60Hz 30VA  | AC100V 50/60Hz 30VA  |
|                  | Temperature            | Operating temperature range: 10 to 40 degrees Celsius  | Operating temperature range: 10 to 40 degrees Celsius  |
|                  |                        | Storage temperature range: -10 to 60 degrees Celsius   | Storage temperature range: -10 to 60 degrees Celsius   |
|                  | Unit weight            | 3.1 kg   | 3.1 kg   |
|                  | Dimensions             | 120mm×120mm×240mm  | 120mm×120mm×240mm  |
|                  | Accessories            | 1.AC power adapter set 2.Probe and testing dish 3.Operation manual (CD-ROM)                            |  |
|                  | Options                | 1.Measurement and recording software NEVA CALL 2.Communication cord 3. Ø2mm, Ø 1mm ,Ø 0.5mm flat prove |  |

<sup>\*</sup>Above specifications are subject to changes without notifications.

### Manufactured and distributed by



9-8, Norimatsuhigashi 1-chome, Yahatanishi-ku, Kitakyushu-city, Fukuoka, 807-0837, Japan TEL +81-93-691-4466 FAX +81-93-602-3065

E-mail : info@iiw.co.jp URL : http://iiw-future.com/